

Capillary Counterdiffusion Kits - Data Sheet - SODIUM CHLORIDE pH 4-9

| Precipitant pH 4 | Capillary / Diammeter | Protein Solution | Observations and Pictures |
|---|------------------------------|-------------------------|----------------------------------|
| Sodium Acetate 0.1M Sodium Chloride 3M | # 1 Ø = | | |
| | # 2 Ø = | | |
| | # 3 Ø = | | |
| | # 4 Ø = | | |
| Precipitant pH 5 | Capillary / Diammeter | Protein Solution | Observations and Pictures |
| Sodium Acetate 0.1M Sodium Chloride 3M | # 1 Ø = | | |
| | # 2 Ø = | | |
| | # 3 Ø = | | |
| | # 4 Ø = | | |
| Precipitant pH 6 | Capillary / Diammeter | Protein Solution | Observations and Pictures |
| Tris-HCl 0.1M Sodium Chloride 3M | # 1 Ø = | | |
| | # 2 Ø = | | |
| | # 3 Ø = | | |
| | # 4 Ø = | | |
| Precipitant pH 7 | Capillary / Diammeter | Protein Solution | Observations and Pictures |
| Tris-HCl 0.1M Sodium Chloride 3M | # 1 Ø = | | |
| | # 2 Ø = | | |
| | # 3 Ø = | | |
| | # 4 Ø = | | |
| Precipitant pH 8 | Capillary / Diammeter | Protein Solution | Observations and Pictures |
| Tris-HCl 0.1M Sodium Chloride 3M | # 1 Ø = | | |
| | # 2 Ø = | | |
| | # 3 Ø = | | |
| | # 4 Ø = | | |
| Precipitant pH 9 | Capillary / Diammeter | Protein Solution | Observations and Pictures |
| Tris-HCl 0.1M Sodium Chloride 3M | # 1 Ø = | | |
| | # 2 Ø = | | |
| | # 3 Ø = | | |
| | # 4 Ø = | | |