



Triana Science & Technology S.L.

Parque Tecnológico Ciencias de la Salud, Avenida Innovación, 1
18100 Armilla - Granada (Spain)

Tel: +34 958 750 583 - Fax: +34 958 750 584

Email: [triana@trianatech.com](mailto: triana@trianatech.com) WebSite: <http://www.trianatech.com>

CRYSTALLIZATION SCREEN FORMULATION

SODIUM FORMATE+MgCl ₂	Precipitant 1	[Precipitant 1]	Precipitant 2	[Precipitant 2]	[Buffer]	pH
Box- 1	Sodium Formate	5 M	MgCl ₂	0.4 M	100 mM	4.00
Box- 2	Sodium Formate	5 M	MgCl ₂	0.4 M	100 mM	5.00
Box- 3	Sodium Formate	5 M	MgCl ₂	0.4 M	100 mM	6.00
Box- 4	Sodium Formate	5 M	MgCl ₂	0.4 M	100 mM	7.00
Box- 5	Sodium Formate	5 M	MgCl ₂	0.4 M	100 mM	8.00
Box- 6	Sodium Formate	5 M	MgCl ₂	0.4 M	100 mM	9.00