



Triana Science & Technology S.L.

Parque Tecnológico Ciencias de la Salud, Avenida Innovación, 1
18100 Armilla - Granada (Spain)
Tel: +34 958 750 583 - Fax: +34 958 750 584
Email: [triana@trianatech.com](mailto: triana@trianatech.com) WebSite: <http://www.trianatech.com>

CRYSTALLIZATION SCREEN FORMULATION

SODIUM MALONATE	Precipitant 1	[Precipitant 1]	Buffer	[Buffer]	pH
Box- 1	Sodium Malonate	3 M	Sodium Acetate	100mM	4.00
Box- 2	Sodium Malonate	3 M	Sodium Acetate	100mM	5.00
Box- 3	Sodium Malonate	3 M	Tris-HCl	100 mM	6.00
Box- 4	Sodium Malonate	3 M	Tris-HCl	100 mM	7.00
Box- 5	Sodium Malonate	3 M	Tris-HCl	100 mM	8.00
Box-6	Sodium Malonate	3 M	Tris-HCl	100 mM	9.00